




How to Handle Your Ocular Prosthesis

Removing the Prosthesis



Pull down on the lower lid with your index finger. You may be able to place your finger under bottom of the prosthesis.


Look up, and allow the prosthesis to slide out over the lower lid.

Inserting the Prosthesis


Dot or Arrow signifies the top of the prosthesis.

Pull up on the upper lid to clearly expose the eye socket cavity.

Gently slide the top of the prosthesis under the upper lid.





... with a Suction Cup



Wet the hollow end of the suction cup, squeeze the stem, and press it against the prosthesis.

Relax the squeeze on the stem and confirm that there is suction to the prosthetic eye.

Gently slide the top of the prosthesis under the upper lid.



Hold down the lower lid with your finger as you tilt the prosthesis out and over the lower lid.

Hold the prosthesis in place with your index finger as you pull down the lower lid with your other hand.

